

# Guideline on running Citrullination\_Diagnostic\_Ion\_Analysis

## Requirements

### Python requirements

- Python version 3.6 or greater
- Dependencies: pandas, numpy, pyteomics, itertools, collections, statistics, re, io

### File requirements

- Input search result file: The input file must be a .csv file(s) containing the following three columns (see next page for the details):

'Title'	MS2 spectrum title as written in MGF file
'Peptide'	Peptide sequence with modification delta mass rounded up to third decimal places
'Charge'	Charge state of the peptide
- Spectrum file: The spectrum file must be a Mascot Generic Format (MGF) file(s) containing MS2 spectra corresponding to those matched to the PSMs in the input search result file. If MS2 spectra in the input search file and spectrum file are not equivalent, only the common MS2 spectra will be retained and subsequently processed.

# Input search result file

Example input search result file:

A	B	C
Title	Peptide	Charge
20160312_02_A1.10012.10012.2 File:"20160312_02_A1.raw", NativeID:"controllerType=0 controllerNumber=1 scan=1001 YETSGIGEAR+0.984VK		2
20160312_02_A1.10045.10045.2 File:"20160312_02_A1.raw", NativeID:"controllerType=0 controllerNumber=1 scan=1004 NIVTPR+0.984TPPPSQGK		2
20160312_02_A1.10116.10116.3 File:"20160312_02_A1.raw", NativeID:"controllerType=0 controllerNumber=1 scan=1011 NIVTPR+0.984TPPPSQGK		3
20160312_02_A1.10222.10222.2 File:"20160312_02_A1.raw", NativeID:"controllerType=0 controllerNumber=1 scan=1022 NIVTPR+0.984TPPPSQGK		2
20160312_02_A1.10334.10334.3 File:"20160312_02_A1.raw", NativeID:"controllerType=0 controllerNumber=1 scan=1033 NIVTPR+0.984TPPPSQGK		3
20160312_02_A1.10362.10362.3 File:"20160312_02_A1.raw", NativeID:"controllerType=0 controllerNumber=1 scan=1036 TPSTAHLR+0.984VPK		3
20160312_02_A1.10418.10418.2 File:"20160312_02_A1.raw", NativeID:"controllerType=0 controllerNumber=1 scan=1041 NIVTPR+0.984TPPPSQGK		2
20160312_02_A1.10479.10479.2 File:"20160312_02_A1.raw", NativeID:"controllerType=0 controllerNumber=1 scan=1047 AQSR+0.984EQLAALK		2
20160312_02_A1.1054.1054.2 File:"20160312_02_A1.raw", NativeID:"controllerType=0 controllerNumber=1 scan=1054" DSR+0.984SGSPM+15.995AR		2
20160312_02_A1.10602.10602.2 File:"20160312_02_A1.raw", NativeID:"controllerType=0 controllerNumber=1 scan=1060 NIVTPR+0.984TPPPSQGK		2
20160312_02_A1.10646.10646.2 File:"20160312_02_A1.raw", NativeID:"controllerType=0 controllerNumber=1 scan=1064 Q+0.984KR+0.984LQ+0.984AM+15.995Q+0.984K		2
20160312_02_A1.10671.10671.3 File:"20160312_02_A1.raw", NativeID:"controllerType=0 controllerNumber=1 scan=1067 SGSEAGSPRR+0.984PRRQR		3
20160312_02_A1.1073.1073.2 File:"20160312_02_A1.raw", NativeID:"controllerType=0 controllerNumber=1 scan=1073" R+0.984GGGGRR+0.984SK		2
20160312_02_A1.10764.10764.2 File:"20160312_02_A1.raw", NativeID:"controllerType=0 controllerNumber=1 scan=1076 R+0.984FIN+0.984DMVK		2
20160312_02_A1.10769.10769.2 File:"20160312_02_A1.raw", NativeID:"controllerType=0 controllerNumber=1 scan=1076 NIVTPR+0.984TPPPSQ+0.984GK		2
20160312_02_A1.10874.10874.2 File:"20160312_02_A1.raw", NativeID:"controllerType=0 controllerNumber=1 scan=1087 MAR+0.984EAFFAEQER		2
20160312_02_A1.11026.11026.3 File:"20160312_02_A1.raw", NativeID:"controllerType=0 controllerNumber=1 scan=1102 RGR+0.984PPKDEK		3
20160312_02_A1.11286.11286.2 File:"20160312_02_A1.raw", NativeID:"controllerType=0 controllerNumber=1 scan=1128 N+0.984R+0.984Q+0.984VIC+57.021VTLK		2
20160312_02_A1.11398.11398.2 File:"20160312_02_A1.raw", NativeID:"controllerType=0 controllerNumber=1 scan=1139 GGTSR+0.984ALAAASVK		2
20160312_02_A1.11489.11489.3 File:"20160312_02_A1.raw", NativeID:"controllerType=0 controllerNumber=1 scan=1148 EEFER+0.984Q+0.984N+0.984KQLR		3
20160312_02_A1.11557.11557.2 File:"20160312_02_A1.raw", NativeID:"controllerType=0 controllerNumber=1 scan=1155 TVEMR+0.984DGEVIK		2
20160312_02_A1.11735.11735.2 File:"20160312_02_A1.raw", NativeID:"controllerType=0 controllerNumber=1 scan=1173 Q+0.984Q+0.984IADLR+0.984EDLKR		2


Format for peptide sequences should follow that of MS-GF+ search result. Specifically, modification delta masses should be rounded up to third decimal places. Currently allowed modifications are as follows:

Modification	Mod on peptide
Carbamidomethyl Cys	C+57.021
Oxidation Met	M+15.995
Deamidated Asn	N+0.984
Deamidated Gln	Q+0.984
Citrullinated Arg	R+0.984
Pyro-Glu from Glu	E-17.027
Pyro-Glu from Gln	Q-18.011
iTRAQ 4plex Lys	K+144.102
iTRAQ 8plex Lys	K+304.205
TMT Lys	K+229.163
iTRAQ 4plex N-term	+144.102
iTRAQ 8plex N-term	+304.205
TMT N-term	+229.163
Acetyl N-term	+42.011

## Downloading the Python script

The Python script can be downloaded via following GitHub page:

<https://github.com/Sunghyun-Huh/Citrullination-Diagnostic-Ion-Analysis>




Users can download the Jupyter notebook via the following steps:

- Click on the 'Citrullination\_Diagnostic\_Ion\_Analysis.ipynb'.
- Click on the 'Raw'.
- Press ctrl+s and manually type 'ipynb' after the filename to download as a .ipynb file.

# Running the Python script

A snapshot of the Jupyter notebook:



The screenshot displays a Jupyter Notebook interface. At the top, the title bar reads "jupyter Citrullination\_Diagnostic\_Ion\_Analysis" with a "Last Checkpoint: 20 minutes ago (autosaved)" status. A "Logout" button is in the top right. Below the title bar is a menu bar with "File", "Edit", "View", "Insert", "Cell", "Kernel", "Widgets", and "Help". A toolbar contains icons for file operations, a "Run" button, and a "Code" dropdown menu. The notebook area shows two code cells. The first cell contains import statements for various libraries. The second cell, titled "0. Initial parameters", contains logic for setting fragmentation methods and ion types based on user-defined parameters.

```
In [ ]: import pandas as pd
import numpy as np
import os
import glob
import statistics
import itertools
import re
import pickle
from pyteomics import mgf, mass
from collections import Counter
```

### 0. Initial parameters

```
In [ ]: # Fragmentation method used
Frag_method = 'HCD' # HCD, CID, ETD, ECD

# Ion types for theoretical spectrum generation
ion_type = ['Precun', 'y', 'b', 'a', 'z', 'c', 'INT', 'IM']

# Generate each ion type
if Frag_method == 'HCD' or Frag_method == 'CID':
    annot_yb = True # y-, b-, a-ion
    annot_zc = False # z-, c-ion
elif Frag_method == 'ETD' or Frag_method == 'ECD':
    annot_yb = False # y-, b-, a-ion
    annot_zc = True # z-, c-ion
```

# Running the Python script

A snapshot of initial parameters settings:

## 0. Initial parameters

```
In [ ]: # Fragmentation method used
Frag_method = 'HCD' # HCD, CID, ETD, ECD

# Ion types for theoretical spectrum generation
ion_type = ['Precur', 'y', 'b', 'a', 'z', 'c', 'INT', 'IM']

# Generate each ion type
if Frag_method == 'HCD' or Frag_method == 'CID':
    annot_yb = True # y-, b-, a-ion
    annot_zc = False # z-, c-ion
elif Frag_method == 'ETD' or Frag_method == 'ECD':
    annot_yb = False # y-, b-, a-ion
    annot_zc = True # z-, c-ion
annot_precur = True # precursor ion
annot_INT = True # internal ion
annot_IM = True # immonium ion
annot_dict = {
    'Precur': annot_precur,
    'y': annot_yb,
    'b': annot_yb,
    'a': annot_yb,
    'z': annot_zc,
    'c': annot_zc,
    'INT': annot_INT,
    'IM': annot_IM
}

# MS2 mass tolerance (ppm)
ms2_ppm = 15

# Signal-to-noise (SNR) filter for MS2 spectrum
apply_SNR = True # Apply SNR filter
SNR = 2 # SNR threshold
low = 0.05 # Define low x% intensity as baseline noise level

# Maximum charge state of sequence ions
max_charge = 2 # 2, 3, ... 'max'

# Maximum number of neutral loss from a single ion
max_NL = 3
```

Explanations of initial parameters are as follows:

Frag_method	Fragmentation method used in the input data (value = 'HCD', 'CID', 'ETD', or 'ECD'). If set as 'HCD' or 'CID', y-ion, b-ion, and a-ion will be generated for theoretical spectrum. If set as 'ETD' or 'ECD', z-ion, c-ion will be generated for theoretical spectrum. Commonly, precursor, internal, and immonium ion will be generated for all fragmentation method used.
ms2_ppm	MS2 level mass tolerance in ppm (default = 15 ppm)
apply_SNR	Determine whether to apply signal-to-noise filter to remove noise peaks (value = True or False; default = True)
SNR	Signal threshold level. The average intensity of noise peaks (as defined in 'low') multiplied by this signal threshold level will be the final signal-to-noise filter. If 'apply_SNR' = True, all peaks below the signal-to-noise filter will be removed (default = 2)
low	Proportion of MS2 peaks regarded as noise. If 'apply_SNR' = True, all peaks below this noise level will be treated as noise (default = 0.05)
max_charge	Maximum charge state of fragment ions (default = 2)
max_NL	Maximum number of neutral losses from a single ion (default = 3)

# Running the Python script

A snapshot of codes for loading input files:

## 1. Input files

```
In [ ]: # Set current working directory
PATH = "F:/Project/"
os.chdir(PATH)

In [ ]: # Input files
spec_files = glob.glob('spectrum_file.mgf') # MGF file(s)
search_files = glob.glob('search_result_file.csv') # Search result file(s)
```

Users can upload local input files via the following steps:

- Set the directory in which the input files are located.
- Copy and paste the input filenames. In case of multiple MGF or search result files in the same directory, type in '\*.mgf' or '\*.csv'.

[illegible]

Column	Explanation
mod_peptide	Simplified peptide with a predefined set of symbols for modifications
Pep_length	Peptide length
mz_Precursor	Theoretical precursor $m/z$
Cit_Count	Number of citrullinated sites
Total_NL_label	Annotations of all diagnostic neutral loss ions
precNL_label	Annotations of precursor neutral losses
seqNL_label	Annotations of sequence ion neutral losses
intNL_label	Annotations of internal ion neutral losses
Total_INT_label	Annotations of all diagnostic internal ions
Dipeptide_label	Annotations of diagnostic dipeptides
Tripeptide_label	Annotations of diagnostic tripeptides
IM_NH3_label	Annotation of IM(Cit)-NH <sub>3</sub>
Total_NL_count	Number of all diagnostic neutral loss ions
precNL_count	Number of precursor neutral losses
seqNL_count	Number of sequence ion neutral losses
intNL_count	Number of internal ion neutral losses
Total_INT_count	Number of all diagnostic internal ions
Dipeptide_count	Number of diagnostic dipeptides
Tripeptide_count	Number of diagnostic tripeptides
IM_NH3_count	Number of IM(Cit)-NH <sub>3</sub>
Cit_probability	Probability (P) of citrullination status calculated by the EN model (HCD data only)
Cit_prediction	Classification of citrullination status using a P cutoff >0.5 (HCD data only)